
Product Name: Pro-Lux™ Proctoscope

ZDOCK271 Rev. 09

Language	Page
en	2-5
de	6-9
fr	10-13
nl	14-17
da	18-21
fi	22-25
sv	26-29
no	30-33

ROCKET MEDICAL PLC
Sedling Road, Washington,
NE38 9BZ, UK

ROCKET MEDICAL GmbH
Am Rosengarten 48,
15566 Schöneiche.
Germany

ROCKET MEDICAL GmbH
Am Rosengarten 48,
15566 Schöneiche.
Germany

ROCKET MEDICAL LLC
50 Corporate Park Drive #890.
East Pembroke. MA. 02359.
USA

Australian Sponsor:

ROCKET MEDICAL PTY LTD
Suite 6, 157 Gordon Street,
Port Macquarie.
NSW 2444. Australia a

Rocket, ROCKET, the R device and the Rocket Medical logo are a registered trademark of Rocket Medical Plc.

PRO-LUX™ PROCTOSCOPE INSTRUCTIONS FOR USE (Including Battery Removal)

Product Name: Pro-Lux™ Proctoscope

Product Code: R53905

Device Image:


Image I: R53905 Pro-lux™ Proctoscope

Intended Use: The Rocket Pro-Lux™ Proctoscope is intended to enable visualisation of and provide access to the anus and distal rectum.

Indications: The Rocket Pro-Lux™ Proctoscope is indicated for use for patients who require visualisation of and access to the anus and distal rectum for the diagnosis and treatment of haemorrhoids or for the retrieval of small foreign bodies,

Contraindications:

Absolute Contraindications

- Imperforate Anus
- Immunosuppressed patient
- Absence of anus following surgical excision
- Stricture
- Moderate to severe anal pain
- Prolapsed thrombosed internal haemorrhoids

Relative Contraindications:

- Severe rectal pain may preclude an awake proctoscope procedure and require sedation
- Recent acute myocardial infarction
- Acute abdomen
- Patients with valvular heart disease or prosthetic valves require prophylactic antibiotics before the procedure
- Patients with coagulopathy (or anticoagulation medication) if biopsy is to be considered
- Major rectal trauma
- Post-operative status (anal surgery)

WARNING:	Battery packs must NOT be allowed to come into contact with the batteries of other devices as they may short-out, leak fluid, overheat or explode. No modification of this equipment is allowed. This device contains a battery power source. Do not touch power source and patient simultaneously.
PRECAUTION:	User should wear gloves and any other suitable personal protective equipment (not provided).

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

PRO-LUX™ PROCTOSCOPE

INSTRUCTIONS FOR USE (Including Battery Removal)

Benefit Risk: Rocket Medical PLC concludes that the overall medical benefits of Pro-Lux Proctoscope, outweighs the possible risks when used according to the intended use.

Performance Claims/Characteristics and Clinical Benefit: The intended performance for the proctoscope is for the visualisation and examination of the anus and distal rectum.

Undesirable Side Effects:

- Discomfort post examination
- Tearing of the perianal skin or mucosa
- Abrasion or tearing of hemorrhoidal tissue
- Infection
- Uncomfortable urge to defecate

Description of the device: The Pro-Lux™ Proctoscope is a single use, non-sterile illuminating device for the use of diagnosis and treatment of the anus and rectum. The light-source provides visibility of the area being examined

User: This device should only be used by or under the supervision of trained personnel and in conjunction with current local clinical practice guidelines.

Patient Population: The intended patient demographics are adult patients requiring the examination and evaluation of their anal cavity and distal rectum for diagnosis and treatment of haemorrhoids, and/or retrieval of small foreign bodies.

Environment: This device is for use in a clinical setting.

For Single Use Only: The device may not have structural integrity to support reuse, resulting in harm to the patient and/or user. Re-use may result in harm to patients and users from cross-contamination and infection.

Non-Sterile: The device is supplied non-sterile. They are not intended to be sterilised prior to use. Do not sterilise as doing so may compromise the structural integrity of the device, leading to device failure.

MRI Safety: MRI compatibility has not been verified.

General Information:

Storage: Keep away from sunlight. Keep dry.

ME Classification: Internally Powered ME Equipment.

Applied part Type: Type B Applied Part.

Applied Parts Table

Description	Contact
3.0V lithium manganese button cell CR1616	User - removing battery for disposal
Proctoscope Battery Puller	User - removing battery for disposal
Proctoscope Back Plate	User - holds during procedure and battery removal for disposal
Proctoscope	User - contact during procedure
Battery Pull Tab	User - Remove to activate light source
Proctoscope obturator	User - contact during procedure Patient - contact during procedure
Proctoscope	Patient - contact during procedure

Physical Characteristics:

Mass: 0.047kg

Dimensions: 1300x1700x500mm

Expected battery life: 35 minutes

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

PRO-LUX™ PROCTOSCOPE

INSTRUCTIONS FOR USE (Including Battery Removal)

Operating Conditions:

- Ambient temperature range: +5°C to +40°C
- Relative humidity range: 15% to 93%
- The device must be operated at ambient pressure levels between 70kPa to 106kPa

Suitability for continuous operation: Device (light source) suitable for continuous use up to 30 minutes during a single use procedure

Symbols and safety signs used on the device and labelling:

Symbol or Safety Sign	Meaning	Symbol or Safety Sign	Meaning
	Manufacturer		Keep away from sunlight
	Caution: Federal law restricts this device to sale by or on the order of a physician.		Do not re-use
	Date of manufacture		Consult the instructions for use
	Use-by date		Caution
	Batch code		Does not contain natural rubber latex
	Catalogue number		Do not use if package is damaged and consult instructions for use
	Non-sterile		Medical device
	Distributor		Do not use if package is damaged and consult instructions for use
	CE Mark		UKCA Symbol
	European Representative		Type B Applied Part
	Keep dry		Importer

Operating Instructions:

CHECK: Do not use if the packaging is open or damaged or the device is damaged.

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

PRO-LUX™ PROCTOSCOPE

INSTRUCTIONS FOR USE (Including Battery Removal)

USE:

1. Remove device from its packaging.
2. Remove tab to activate the light source (**Figure 1**).
3. Place the patient in the left lateral decubitus position with knees flexed towards the chest.
4. Inspect anus and surrounding areas for any abnormalities.
5. Apply lubricant to a gloved finger and perform a digital rectum examination.
6. Lubricate the proctoscope with the obturator in place.
7. Spread the buttocks and gently insert the proctoscope into the anal canal. If required, ask the patient to take a few deep, gentle breaths and to bear down slightly to help make the insertion easier. Gently advance the instrument towards the umbilicus until the full length is inserted.
8. After the proctoscope is completely inserted, remove obturator.
9. Remove any faecal matter with a large swab (not supplied).
10. Perform examination, treat haemorrhoids or remove foreign bodies (where applicable).
11. Slowly rotate the proctoscope as it is withdrawn.
12. Once the device is fully withdrawn from the patient's anal canal, wipe away any excess lubricant with tissue.
13. In accordance with local hospital policy, remove batteries before disposal. Push tab into handle to release battery cover (**Figure 2**).
14. Lift cover to expose battery pack ready for removal (**Figure 3**).
15. Pull the red tab shown to remove battery pack (**Figure 4**).
16. To ensure safety, place each battery pack into a separate cell of the disposal tray provided with each carton of Pro-Lux Proctoscopes. Dispose of the filled tray in accordance with local hospital policy.


Figure 1


Figure 2


Figure 3


Figure 4

Disposal: This device should be handled and disposed of in accordance with policy of the healthcare setting and with regard to all applicable regulations, including but without limitation to, those pertaining to human health & safety and care of the environment. Failure to do so may increase the risks of infection or other microbial hazards.

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

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