
Pudendal Block Needle – NRFit®

ZDOCK330 Rev. 02

Language	Page
En	2-6
Da	7-11
De	12-16
Es	17-21
Fr	22-26
It	27-31
Nl	32-36
No	37-41
Sv	42-46

ROCKET MEDICAL PLC
Sedling Road, Washington,
NE38 9BZ, UK

ROCKET MEDICAL GmbH
Am Rosengarten 48,
15566 Schöneiche,
Germany

ROCKET MEDICAL GmbH
Am Rosengarten 48,
15566 Schöneiche,
Germany

ROCKET MEDICAL LLC
50 Corporate Park Drive #890.
East Pembroke. MA. 02359.
USA

ROCKET MEDICAL PTY LTD
Suite 6, 157 Gordon Street,
Port Macquarie.
NSW 2444. Australia

Australian Sponsor:

The R device and the Rocket Medical logo are registered trademarks of Rocket Medical Plc. NRFit® is a trademark of GEDSA and is used with permission.

PUDENDAL BLOCK NEEDLE – NRFit® INSTRUCTIONS FOR USE

Product Name: Pudendal Block Needle – NRFit®

Product Code: R57040-NRFIT

Device Image:


Figure 1 – Labelled Image of Pudendal Block Needle – NRFit® – with Transport Clip in Place

Intended Use: Intended for transvaginal injection of local anaesthetic agents into the pudendal nerve area for the creation of a regional pudendal nerve block during the second stage of labour and in obstetric procedures.

Indications: Indicated for use when a regional pudendal nerve block is required:

- During the second stage of labour,
- Prior to the repair of an episiotomy or perineal laceration, or
- Prior to an instrumental delivery.

Contraindications: Use of the device is contraindicated:

- when fetal head engagement prevents vaginal palpation as the fetus may be injured by the needle or toxicity may arise due to direct injection of anaesthetic into the fetus
- where there is known severe sensitivity to local anaesthetic agents, as this may lead to anaphylaxis
- in the presence of active infection of the skin or injection site, or within the ischiorectal area or the adjacent structures, including the vagina or perineum. This may lead to worsening of the infection.
- in patients with uncorrected coagulation disorders as this may lead to severe bleeding and hematomas.

WARNING: Care must be taken to avoid intravascular administration into the pudendal vessels which lie lateral to the pudendal nerve. Aspirate to confirm the needle placement is not intravascular prior to injection. Appropriate post procedural monitoring of the patient and the fetus is mandatory. Ensure that provision for resuscitation is made, should adverse reactions or toxicity occur. Exercise caution in the combination and dosage of sedatives and narcotics to avoid possible airway obstruction or aspiration. This device is intended for an individual patient, during a single procedure and should be discarded after use. Do not resterilise.

Do not use if packaging is damaged or broken or has been previously opened.

Do not use if the pack label indicates that the device has expired.

Limit use to trained, experienced or supervised operators.

CAUTION: Ensure transport clip is still in place prior to repositioning of finger support ring.

To minimise the risk of needlestick injury, hold the connector when removing the transport clip to ensure that the needle does not protrude from the needle guide.

When guiding the needle to the injection site, take care to ensure that the needle is sheathed in the needle guide.

Discard the removable finger support ring and place the needle and sleeve assembly directly into an approved sharps container.

Do NOT attempt to reattach the transport clip prior to disposal.

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pnrf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

Benefit Risk: Rocket Medical plc has taken all necessary steps to ensure that residual risks associated with the use of Pudendal Block Needle are reduced as far as possible through application of existing state of the art techniques in the design and manufacture of these medical devices to ensure safe usage. Rocket Medical plc concludes that the overall medical benefits of Pudendal Block Needle outweigh the possible risks when used according to the intended use.

Performance Claims/Characteristics and Clinical Benefit:

Performance Claim: To provide transvaginal access to the pudendal nerve space for delivery of local anaesthetic to create a regional pudendal nerve block.

Clinical Benefit (Indirect): Rocket Medical Pudendal Block Needle facilitates the administration of anaesthetics and does not have a direct therapeutic or diagnostic function. The clinical benefit is indirect (MDCG 2020-6).

Undesirable Side Effects:

- Pain at the injection site.
- Laceration of the vaginal mucosa due to accidental injury from the needle
- Hematoma from perforation of a blood vessel during needle insertion.
- Infection along the psoas muscle or obturator internus.
- Ischial region paraesthesia resulting from nerve injury.
- Sacral neuropathy resulting from nerve injury.
- Systemic toxicity due to intravascular injection.

Description of the device: The pudendal block needle consists of a 20G x 130mm needle, needle guide and a connector hub to allow a prefilled NRFit® lock syringe to be attached. The needle guide prevents over insertion of the needle into the tissue during use, facilitates access and offers the ability to withdraw the needle to prevent unintended needle injuries. The support ring can be placed in two different positions depending on the approach and clinical preference. The device also has a clip to stop the needle protruding during transport before use.

The device is inserted vaginally for the introduction of a local anaesthetic agent to create a regional pudendal nerve block for analgesia. The depth of submucosal penetration is limited by the device to minimise the risk of injury to the vagina and the fetus.

User: The intended user is a healthcare professional proficient in the transvaginal injection of anaesthetic agents to create a regional pudendal nerve block, or a person in training under supervision of a suitable competent healthcare professional, working in accordance with local and national guidelines.

Patient Population: The intended patient demographic is patients requiring a transvaginal pudendal nerve block for analgesia during the second stage of labour and in obstetric procedures.

Environment: Clinical setting, with access to sharps bin.

For Single Use Only: The device may not have structural integrity to support reuse, resulting in harm to the patient and/or user. Re-use may result in harm to patients and users from cross-contamination and infection.

Sterile: The device(s) is supplied sterile by Ethylene Oxide (EO). Do not re-sterilise; it may not be possible to achieve the required sterility assurance level upon re-sterilisation of a used device. Re-sterilisation may also compromise the structural integrity of the device, leading to device failure.


MRI Safety: This device is not intended to be routinely found in an MRI environment and compatibility has not been verified.

General Information:

- The device is compatible with syringes with an NRFit® connector and is NOT compatible with Luer lock syringes.

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

Symbols and safety signs used on the device and labelling:

Symbol or Safety Sign	Meaning	Symbol or Safety Sign	Meaning
	Manufacturer		Keep away from sunlight
	Authorized representative in the European Community/European Union		Keep dry
	Caution: Federal law restricts this device to sale by or on the order of a physician.		Do not re-use
	Date of manufacture		Consult the instructions for use
	Use-by date		Does not contain natural rubber latex
	Batch code		Medical device
	Catalogue number		Do not use if package is damaged and consult instructions for use
	Sterilized using ethylene oxide		Single sterile barrier system
	Do not re-sterilise		Importer
	Single sterile barrier system with protective packaging outside		Distributor

Operating Instructions:

CHECK:

- Do not use if the packaging is open or damaged or the device is damaged.
- Ensure that you have an NRFit® syringe, NRFit® drawing up needle and the required anaesthetic.
- Ensure that suitable equipment is available for post-procedural maternal and fetal monitoring.
- Provision should be made for resuscitation in the event that there is an adverse reaction to the anaesthetic or should toxicity occur.

PREPARE:

Using aseptic technique, remove the device from the outer bag.

The finger support ring can be located in one of two positions (see figure 1) depending on the approach and surgical preference.

CAUTION: Ensure transport clip is still in place prior to repositioning of finger support ring.

Following repositioning (if required), remove the transport clip from the needle.

CAUTION: To minimise the risk of needlestick injury, hold the connector when removing the transport clip to ensure that the needle does not protrude from the needle guide.

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.


Figure 2 – Image of Pudendal Block Needle – NRFit® – with Transport Clip Partially Removed

CONTROLS: The depth of submucosal penetration is limited to 12mm by the device to minimise the risk of injury to the vagina and the fetus.

USE:

1. Place the patient in the lithotomy position. Following local hospital procedure, prep the perineum and vagina using an approved antibacterial solution.
2. The ischial spines can be palpated at 4 and 8 o'clock about a finger's length into the vagina. The pudendal nerve runs behind the sacrospinous ligament, about 1cm antero-medial and postero-medial to the ischial spine.
3. Draw up the required anaesthetic into an NRFit® syringe. Attach the NRFit® syringe of local anaesthetic to the needle connector, taking care to ensure that the needle is sheathed in the needle guide.
4. To perform a left-sided block, palpate the ischial spine with the index finger of the left hand. Holding the syringe in the right hand, guide the needle/needle guide between the index and middle finger of the left hand toward the ischial spine, taking care to ensure that the needle is sheathed in the needle guide.
5. Triple Injection Technique:
 - a. Place the end of the guide beneath the tip of the ischial spine and advance the needle into the vaginal mucosa.
 - b. Aspirate to ensure that the injection is not intravascular.
 - c. Raise a mucosal wheal with local anaesthetic.
 - d. Advance the needle through the vaginal mucosa until it touches the sacrospinous ligament, which is located approximately 1 cm medial and posterior to the ischial spine.
 - e. Infiltrate the tissue with local anaesthetic.
 - f. Advance the needle further through the sacrospinous ligament for a distance of 1 cm until a loss of resistance is detected. The needle tip now lies in the area of the pudendal nerve.
WARNING: Care must be taken to avoid intravascular administration into the pudendal vessels which lie lateral to the pudendal nerve. Aspirate to confirm the needle placement is not intravascular prior to injection.
 - g. Inject local anaesthetic solution into this region.
 - h. Withdraw the needle into the guide and move the tip of the guide to just above the ischial spine. Reinsert the needle through the mucosa.
 - i. Aspirate to ensure that the injection is not intravascular.
 - j. Infiltrate the tissue with local anaesthetic.
 - k. Sheath the needle by pulling it into the needle guide. Withdraw the needle and needle guide together.
6. Single Injection Technique:
 - a. An alternative technique to the triple injection technique uses a single injection of a larger volume of local anaesthetic when the needle is introduced nearly 1cm through the sacrospinous ligament, medial and posterior to the ischial spine.
 - b. **WARNING: Care must be taken to avoid intravascular administration into the pudendal vessels which lie lateral to the pudendal nerve. Aspirate to confirm the needle placement is not intravascular prior to injection.**
 - c. Infiltrate the tissue with local anaesthetic.
 - d. Sheath the needle by pulling it into the needle guide. Withdraw the needle and needle guide together.
7. To perform a right-sided block, repeat these steps holding the syringe in the left hand and guiding the needle/needle guide between the index and middle finger of the right hand toward the ischial spine taking care to ensure that the needle is sheathed in the needle guide.
8. The patient must be checked bilaterally to establish the extent and efficacy of the block before proceeding with the surgical procedure. The block typically becomes effective 5-20 minutes after infiltration and may last for 20-60 minutes depending on the agent used.

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

9. Appropriate post procedural monitoring of the patient and the fetus is mandatory. Ensure that provision for resuscitation is made, should adverse reactions or toxicity occur. Exercise caution in the combination and dosage of sedatives and narcotics to avoid possible airway obstruction or aspiration.

When an early pudendal block is given, the block may later be repeated as long as the maximum dose of local anaesthetic is not exceeded.

Adapted from Hemant K. Satpathy: *Transvaginal Pudendal Nerve Block* 12/2020; <https://emedicine.medscape.com/article/83078-overview>


Disposal: This device and the consumables used with it, should be handled and disposed of in accordance with policy of the healthcare setting and with regard to all applicable regulations, including but without limitation to, those pertaining to human health and safety and care of the environment. Failure to do so may increase the risks of infection or other microbial hazards. Take care when handling sharps, to avoid needlestick injuries. Discard the removable finger support ring and place the needle and sleeve assembly directly into an approved sharps container. Do not attempt to reattach the transport clip prior to disposal.

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

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