
Product Name: Embryo Transfer Catheters

ZDOCK067 Rev. 28

Language	Page
en	2-5
es	6-9
de	10-13

ROCKET MEDICAL PLC
Sedling Road, Washington,
NE38 9BZ, UK

Australian Sponsor:

ROCKET MEDICAL GmbH
Am Rosengarten 48,
15566 Schöneiche,
Germany

ROCKET MEDICAL GmbH
Am Rosengarten 48,
15566 Schöneiche,
Germany

ROCKET MEDICAL LLC
50 Corporate Park Drive #890,
East Pembroke. MA. 02359,
USA

ROCKET MEDICAL PTY LTD
Suite 6, 157 Gordon Street,
Port Macquarie,
NSW 2444. Australia

Rocket, ROCKET, EMBRYON, the R device and the Rocket Medical logo are registered trademarks of Rocket Medical Plc.

EMBRYON® BULBTIP™ EMBRYO TRANSFER CATHETERS

INSTRUCTIONS FOR USE

Product Name: Embryon® BulbTip™ Embryo Transfer Catheters and Embryo Transfer Trial Catheters.

Product Codes: R57635-00-18, R57635-00-23, R57635-EC-18, R57635-EC-23 R57636-00-18 and R57636-00-23.

Device Image:

Figure 1: (A) Embryon® BulbTip™ Embryo Transfer Catheter & (B) Embryon® BulbTip™ Echogenic Embryo Transfer Catheter

Intended Use:

Embryo Transfer Catheters- For the introduction of human embryos into the uterine cavity during assisted reproduction.

Trial Catheters - Embryo Trial Transfer Catheters are intended for the trial procedure prior to the introduction of human embryos into the uterine cavity during assisted reproduction.

Indications:

Embryo Transfer Catheters - Rocket Medical Embryo Transfer Catheters are indicated for the placement of human embryos into the uterine cavity for patients requiring assisted reproduction.

Trial Catheters - Rocket Medical Embryo Trial Transfer Catheters are indicated for a trial procedure to ensure correct placement prior to the transfer of the human embryo into the uterine cavity for patients requiring assisted reproduction.

Contraindications:

The devices are contraindicated for:

- **Infection:** if embryo transfer is performed local infection can be spread into the peritoneal cavity causing systemic infection pelvic infection and risk to future fertility.
- **The presence of, or after recent Pelvic Inflammatory Disease:** if embryo transfer is performed infection can be spread causing systemic infection.
- **Intra-fallopian procedures:** embryo transfer catheters are not designed to enter the fallopian tubes and could cause tubal damage.

WARNINGS:	Cervical Assessments should be made prior to procedure Not for long term use. Please be aware that withdrawal of the inner catheter through the outer sheath while the tip of the catheter is in the uterus or endocervical canal may result in detachment of the echogenic band.
------------------	---

Precaution:	Device is intended to be stored and used in environments where the temperature falls between 5°C - 30°C
--------------------	--

Transfer Catheters are reduced as far as possible through application of existing state of the art techniques in the design and manufacture of these medical devices to ensure safe usage. Rocket Medical plc concludes that the overall medical benefits of Embryo Transfer Catheters, outweigh the possible risks when used according to the intended use.

Performance Claims and Clinical Benefit: Embryo Transfer Catheters are used to place a human embryo into the uterine cavity as part of the ART procedure.

Undesirable Side Effects: Infection, pain and a risk of multiple births may occur if the patient has been on fertility hormone treatment, pain and infection. Accidental detachment of the echogenic band from the inner has been reported in a small percentage of cases.

Description of the device: The device is constructed of two primary components

1. Soft inner catheter
2. Outer sheath with 5 x 1cm markings from the tip to indicate depth of penetration into the endocervix

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

EMBRYON® BULBTIP™ EMBRYO TRANSFER CATHETERS

INSTRUCTIONS FOR USE

The transparent finger grip has 5 ribs moulded at 1cm intervals. When the inner catheter is inserted into the outer sheath, the coloured marker bead will indicate the position of the catheter tip in relation to the end of the outer sheath. Position of the coloured bead corresponds with:

- a) Proximal edge - tip protrusion = 0cm (flush with the end of sheath)
- b) 1st proximal rib - tip protrusion = 1cm
- c) 2nd proximal rib - tip protrusion = 2cm.
- d) Catheter fully inserted - tip protrusion = 4cm

BulbTip™ Echogenic Transfer catheters have an echogenic band applied 2mm from the distal tip. This echogenic marker may aid placement with ultrasound.

User: These devices may only be used by, or under the supervision of, appropriately trained personnel in embryo transfer procedure, in conjunction with clinical practice guidelines and regulations.

Patient Population: The intended patient demographics is anybody suffering from infertility and wishes to have a child using these methods as an alternative to natural conception.


Environment: Clinical setting

For Single Use Only: the device may not have structural integrity to support reuse, resulting in harm to the patient and/or user. Re-use may result in harm to patients and users from cross-contamination and infection.

Sterile: the devices are supplied sterile by Ethylene Oxide (EO). Do not re-sterilise; it may not be possible to achieve the required sterility assurance level upon re-sterilisation of a used device. Re-sterilisation may also compromise the structural integrity of the device, leading to device failure.

MRI Safety: This device is not intended to be routinely found in an MRI environment and compatibility has not been verified. However, the only metal used in the construction of the device is present to facilitate imaging and the type and quantity is well evidenced to be safe to use with MRI.

Symbols and safety signs used on the device and labelling:

Symbol or Safety Sign	Meaning	Symbol or Safety Sign	Meaning
	Manufacturer		Keep away from sunlight
	Authorized representative in the European Community/European Union		Keep dry
	Temperature limit between 5°C - 30°C		Do not re-use
	Caution: Federal law restricts this device to sale by or on the order of a physician.		Consult the instructions for use
	Date of manufacture		Caution
	Use-by date		Does not contain natural rubber latex

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

EMBRYON® BULBTIP™ EMBRYO TRANSFER CATHETERS

INSTRUCTIONS FOR USE

	Batch code		Medical device
	Catalogue number		Do not use if package is damaged and consult instructions for use
	Sterilized using ethylene oxide		Single sterile barrier system
	Single sterile barrier system with protective packaging outside		Do not resterilise
	Distributor		Importer

Operating Instructions:

CHECK: do not use if the packaging is open or damaged or the device is damaged.

Cervical Assessment: **Rocket Medical Plc** strongly recommends that the position and patency of the cervix and canal is assessed 2-3 weeks before the date of initial transfer. If more than 3 months has elapsed since the previous transfer or assessment, then a further examination should be considered. The total uterine length, including endocervix, should be established by uterine sounding or ultrasonography at the initial patient consultation.

Place the patient in the lithotomy position and expose the cervix using a conventional Cusco or Sims speculum. The procedure can normally be carried out without medication, although light sedation may be necessary in very anxious patients.

Assess the position of the cervix and establish the position of the external cervical os. Swab away any excess mucus with normal saline solution. Using aseptic technique, prepare a matching **Trial Transfer Catheter** and follow one of the two techniques below:

- No significant anteversion or retroversion, cervical os fully visible:** Assemble the inner trial catheter into the outer sheath and secure the luer fitting into the clear outer handle. Locate the distal tip of the trial catheter into the external cervical os and with a very gentle rotating motion, advance the catheter into the uterus. The patient may report slight discomfort when the catheter passes through the internal os and into the endometrial cavity. If the trial transfer catheter will not pass freely into and through the endocervical canal, follow the technique below.
- Anteversion, retroversion with cervical os not fully visible or previous catheter insertion difficult or obstructed:** Remove the inner trial catheter from the outer sheath. Gently grip the outer sheath between gloved forefinger and thumb and *gently* increase the curve to match the amount of anteversion or retroversion. Do not bend sharply as this will kink the outer sheath.

A matching stylet (R57591-ET-18 or R57591-ET-23) may also be inserted into the outer sheath to increase rigidity and allow better access to the severely retro and anteverted uterus. Identify the position of cervical os and advance the outer sheath 2 x 1cm marking. Insert the inner catheter into the clear handle and advance until the inner catheter passes up to the internal os.

Observing the relative position of the marker bead and the moulded finger grooves will indicate the catheter tip position, in 4c, from the distal end of the outer sheath.

Patients who have a fully obstructed cervix or impassable endocervical canal should be considered for cervical surgery or alternative fertility treatment before further attempts at endocervical transfer are made.

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

EMBRYON® BULBTIP™ EMBRYO TRANSFER CATHETERS

INSTRUCTIONS FOR USE

Procedure:

Drawing up embryos: Select a transfer catheter to suit the needs of the patient and surgical preference.

- a) Using aseptic technique, remove the device from the blister tray or peel pouch and attach a sterile microsyringe. Do not remove the inner catheter from the outer sheath as this provides additional support and control.
- b) Under microscopic control, draw 100µl of prepared medium to fill the catheter. Draw up a further 30µl containing the embryos for transfer. The embryos should be retained in the distal third of the catheter until the moment of injection into the uterus. Remove the microsyringe, expel any residual air and fill with medium ensuring that no air remains in the syringe.

OR:

- c) Under microscopic control, draw 50µl of prepared medium into the distal third of the catheter. Draw up a 10-20µl of air. Draw up a further 30µl containing the embryos for transfer. The embryos should be retained in the distal third of the catheter until the moment of injection into the uterus. Remove the microsyringe and fill with medium ensuring that no air remains in the syringe.
- d) Maintaining aseptic technique, transfer the complete assembly to the operating room immediately prior to transfer.

Transferring the embryos:

- a) Ensure that trial transfer and cervical assessment have been completed before attempting a 'live' transfer. Place the patient in the lithotomy position and expose the cervix using a conventional Cusco or Sims speculum. The procedure can normally be carried out without medication, although light sedation may be necessary in very anxious patients.
- b) Transvaginal or transabdominal ultrasound may be useful in establishing correct placement of the catheter tip in the uterus.
- c) Assess the position of the cervix and establish the position of the external cervical os. Swab away any excess mucus with normal saline solution or prepared medium.
- d) Follow the techniques described for the placement of the Trial Transfer Catheter to place the loaded embryo transfer catheter.
- e) Once in position, inject approximately 20µl of medium to expel the embryos.
- f) Whilst maintaining pressure on the syringe plunger, gently withdraw the catheter from the uterus and return it to the laboratory to check microscopically that all the embryos have been expelled.
- g) The patient may then be returned to the ward for subsequent discharge.

Disposal: This device, its accessories and the consumables used with it, should be handled and disposed of in accordance with policy of the healthcare setting and with regard to all applicable regulations, including but without limitation to, those pertaining to human health and safety and care of the environment. Failure to do so may increase the risks of infection or other microbial hazards. Take care when handling sharps, to avoid needlestick injuries. Ensure sharps are disposed of in sharps containers.

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

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