Product Names: Seldinger Chest Drain Kit

ZDOCKI49 Rev. 22

Language	Page
en	2-7
cn	8-12
hr	13-17
da	18-22
nl	23-28
fi	29-33
fr	34-39
de	40-45
el	46-51
hu	52-57
it	58-63
no	64-69
pl	70-75
sl	76-80
es	81-86
sv	87-92
tr	93-98
vi	99-104

ROCKET MEDICAL PLC
Sedling Road, Washington,
NE38 9BZ, UK

ROCKET MEDICAL GmbH
Am Rosengarten 48,
15566 Schöneiche.
Germany

ROCKET MEDICAL GmbH
Am Rosengarten 48,
15566 Schöneiche.
Germany

ROCKET MEDICAL LLC
50 Corporate Park Drive #890.
East Pembroke. MA. 02359.
USA

Australian Sponsor:

Rocket Medical Pty Ltd
Suite 6, 157 Gordon Street,
Port Macquarie.
NSW 2444. Australia

Rocket, ROCKET, the R device and the Rocket Medical logo are registered trademarks of Rocket Medical Plc.

Seldinger Chest Drain Kit
INSTRUCTIONS FOR USE

Product Name: Seldinger Chest Drain Kit

Procedure Pack Product Codes: R54544-12-00 and R54544-18-00

Device Image:


Image 1: R54544-12-00


Image 2: R54544-18-00

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

Device Contents: Check you have the full contents of the kit prior to use.

DEVICE REF		DEVICE NAME	QTY	LEGAL MANUFACTURER DETAILS
R54544-12-00	ZPACK016	12Fg Seldinger Catheter	1	ROCKET MEDICAL PLC Sedling Road, Washington, NE38 9BZ, UK www.rocketmedical.com 
		14Fg Dilator with Depth Ring	1	
		Seldinger Widebore Connector	1	
		Cap Sealed – use on Stopcock	1	
		Stepped Male Luer Lock	1	
		Echo Seldinger Needle	1	
		1.95mm Dia LDPE Support Rod	1	
		4 Way Stopcock Large Bore	1	
		Guidewire 70cm Stainless Steel	1	
	ZBLAD013	No 11 Blade Scalpel	1	SWANN MORTON LIMITED , Owlerton Green, Sheffield, S6 2BJ, UK
	ZSYRG015	10ml Luer Lock Syringe	1	BECTON DICKINSON S.A. Crta. Mequinenza s/n, 22520 Fraga, Huesca, España
R54544-18-00	ZPACK019	14Fg Dilator with Depth Ring	1	ROCKET MEDICAL PLC Sedling Road, Washington, NE38 9BZ, UK www.rocketmedical.com 
		Seldinger Widebore Connector	1	
		16Fg Dilator with Depth Ring	1	
		20Fg Dilator with Depth Ring	1	
		18Fg Seldinger Catheter	1	
		Cap Sealed – use on Stopcock	1	
		Stepped Male Luer Lock	1	
		Echo Seldinger Needle 18g	1	
		4 Way Stopcock Large Bore	1	
	ZBLAD013	No 11 Blade Scalpel	1	SWANN MORTON LIMITED , Owlerton Green, Sheffield, S6 2BJ, UK
	ZSYRG015	10ml Luer Lock Syringe	1	BECTON DICKINSON S.A. Crta. Mequinenza s/n, 22520 Fraga, Huesca, España

Intended Use: The products are intended for use by physicians trained and experienced in diagnostic and interventional techniques. It is recommended that British Thoracic Society Guidelines for the insertion of chest drainage tubes are followed.

Indications: The product has been designed for percutaneous introduction of a chest tube for drainage of air and/or pleural fluid.

Contraindications:

- For drainage of pleural space only
- 12fg – not for drainage of thick fluids.

WARNING:	- Over insertion of the dilator into the chest cavity is associated with serious injury and MUST be avoided. - Only insert the dilator sufficient to dilate the superficial track. - This device should only be used by, or under the supervision of, appropriately trained personnel and in conjunction with current local clinical practice guidelines
Caution	Insertion should be just ABOVE the rib to avoid damage to the intercostal neurovascular bundle.
Caution	The wire guide should pass through the needle and advance into the pleural space without resistance.
Caution	Ensure sufficient length of wire guide remains external to the chest to facilitate controlled introduction of the dilators and drain.

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

Caution

- It is important to advance the chest tube assembly into the pleural space whilst maintaining its direction in the same plane as the wire guide. This will make introduction easier and avoid kinking of the wire guide.

Caution

Precaution:

- Ensure all side holes of the chest tube are positioned within the pleural space.
- Potential for DEHP Exposure from use of this device.** DEHP may leach from this device, which may cause reproductive harm in male neonates, pregnant women carrying male fetuses, and peripubertal males. The extent of exposure largely depends upon the medical treatments administered and the duration of the treatment. Lack of research in humans means it is difficult to predict the adverse effects of DEHP because certain animal models may not apply to humans. Most studies utilise mice and rats, however breakdown of DEHP in the human body differs to the mechanism observed in these animals. These effects are also seen only at levels far in excess of normal human exposure. (For further details please see 'POL58' at www.rocketmedical.com.)

Benefit Risk: Rocket Medical PLC concludes that the overall medical benefits of Seldinger Chest Drain Kit, outweigh the possible risks when used according to the intended use.

Performance Claims and Clinical Benefit: The Rocket Seldinger Chest Drainage Set is indicated for the percutaneous insertion of a pleural drain for fluid and air. The Seldinger range is more comfortable for the patient utilising smaller gauge catheters than similar thoracic types.

Undesirable Side Effects: Infection, Pain, Subcutaneous emphysema.

User: The intended user is a healthcare professional proficient in the undertaking of percutaneous procedures, working in accordance with local and national guidelines.

Patient Population: The intended patient demographics are adult patients requiring drainage of air and pleural fluid from the pleural cavity.

Environment: Clinical setting with sterile field with close access to sharps bin e.g. within 2 metres.

For Single Use Only: the device may not have structural integrity to support reuse, resulting in harm to the patient and/or user. Re-use may result in harm to patients and users from cross-contamination and infection.

Sterile: the devices are supplied sterile by Ethylene Oxide (EO). Do not re-sterilise; it may not be possible to achieve the required sterility assurance level upon re-sterilisation of a used device. Re-sterilisation may also compromise the structural integrity of the device, leading to device failure.

MRI Safety: This device is not intended to be routinely found in an MRI environment and compatibility has not been verified. However, the radiopaque material used in the construction of the device is present to facilitate imaging and the type and quantity is well evidenced to be safe to use with MRI.


General Information:

- Devices are supplied with a series of luer lock connectors for attachment to chest drain bottles.

Symbols and safety signs used on the device and labelling:

Symbol or Safety Sign	Meaning	Symbol or Safety Sign	Meaning
	Manufacturer		Keep away from sunlight
	Caution: Federal law restricts this device to sale by or on the order of a physician.		Keep dry
	Date of manufacture		Do not re-use

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

	Use-by date		Consult the instructions for use
	Batch code		Caution
	Catalogue number		Does not contain natural rubber latex
	Sterilized using ethylene oxide		Contains or presence of phthalate
	Do not resterilise		Medical device
	Authorized representative in the European Community/European Union		Do not use if package is damaged and consult instructions for use
	Single sterile barrier system with protective packaging outside		Single sterile barrier system
	Distributor		Importer


Operating Instructions:

CHECK: do not use if the packaging is open or damaged or the device is damaged. Ensure you have the full contents of the kit prior to use (see 'Device Contents').

USE:

1. Following local hospital policy, prepare the chest tube insertion site with an approved solution such as povidone-iodine and drape to maintain and aseptic technique.
2. Identify and mark the insertion point following BTS Guidelines, with specific attention paid to the 'safe-triangle'.
3. After injection of local anaesthetic, make a small skin incision 4-5mm.

CAUTION: Insertion should be just ABOVE the rib to avoid damage to the intercostals neurovascular bundle.

The introducer needle has changed from a Tuohy needle and is now a straight Seldinger needle as below. When removing the protective sheath, hold the silicone ring to prevent it from being removed from the needle.

The measurement provided by using the silicone ring can be used to set the position of the sliding stop on each dilator to suit the size of the patient and help prevent over insertion of the dilators.

4. Attach the introducer needle to a syringe filled with a local anaesthetic e.g. lignocaine, and administer additional local anaesthetic to the skin and underlying tissues.
5. Advance the introducer needle over the *superior* border of the rib and into the pleural space. Fluid or air should be aspirated to verify correct position. The needle should be introduced and directed with appropriate orientation inferiorly or superiorly.
6. Remove the syringe and advance the soft “J” end of the wire guide through the needle and into the pleural space. Before removing the needle, move the silicone measurement ring down to the patient’s skin.

CAUTION: The wire guide should pass through the needle and advance into the pleural space without resistance.

7. Remove the introducer needle.
8. Using the measurement obtained using the depth marker on the introducer needle, set the position of the sliding stop on each dilator to suit the size of the patient.

WARNING:


- Over insertion of the dilator into the chest cavity is associated with serious injury and **MUST** be avoided.
- Only insert the dilator sufficient to dilate the superficial track.

CAUTION: Ensure sufficient length of wire guide remains external to the chest to facilitate controlled introduction of the dilators and drain.

9. While maintaining the wire guide position, dilate the track and opening into the pleural space by advancing and withdrawing, one by one, the supplied dilator(s) over the wire guide. Introduction into the pleural space is facilitated by rotating and advancing the dilators in the same plane as the wire guide to prevent kinking.
10. Remove the final dilator from the wire guide
11. Locate the drain set. Fully insert the stiffening tube, tighten the valve if required to ensure the stiffening tube remains fully inserted to the distal plug. Insert the free end of the guidewire into the distal catheter plug ensuring it travels through the stiffening tube.
12. With the wire guide still positioned within the pleural space, advance the chest tube inserter/chest tube assembly over the wire guide and into the pleural space.


CAUTION:

- It is important to advance the chest tube assembly into the pleural space whilst maintaining its direction in the same plane as the wire guide.
- This will make introduction easier and avoid kinking of the wire guide.

CAUTION: Ensure all side holes of the chest tube are positioned within the pleural space.

- I3. Remove the wire guide and chest tube inserter leaving the chest tube in place.
- I4. The chest tube can now be sutured to the skin and is ready for connection to a suitable sealed drainage system such as the Rocket® Blue Bottle or Rocket® Ambulatory Bag
- I5. Confirm the correct positioning of the drain with imaging as soon as possible.

Disposal: This device, its accessories and the consumables used with it, should be handled and disposed of in accordance with policy of the healthcare setting and with regard to all applicable regulations, including but without limitation to, those pertaining to human health and safety and care of the environment. Failure to do so may increase the risks of infection or other microbial hazards. Take care when handling sharps, to avoid needlestick injuries. Ensure sharps are disposed of in sharps containers.

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