
Product Name: Rocket IPC Replacement Valve

ZDOCKI90 Rev. 14

Language	Page
en	2-4
ar	5-7
cn	8-10
cz	11-13
da	14-16
de	17-19
el	20-23
es	24-26
fi	27-29
fr	30-32
he	33-35
it	36-38
lt	39-41
nl	42-44
pl	45-47
pt	48-50
sl	51-53
sv	54-56
tr	57-59

ROCKET MEDICAL PLC
Sedling Road, Washington,
NE38 9BZ, UK

ROCKET MEDICAL GmbH
Am Rosengarten 48,
15566 Schöneiche.
Germany

ROCKET MEDICAL GmbH
Am Rosengarten 48,
15566 Schöneiche.
Germany

ROCKET MEDICAL LLC
50 Corporate Park Drive #890.
East Pembroke. MA. 02359.
USA

Australian Sponsor:

ROCKET MEDICAL PTY LTD
Suite 6, 157 Gordon Street,
Port Macquarie.
NSW 2444. Australia

Rocket, ROCKET, the R device and the Rocket Medical logo are a registered trademark of Rocket Medical Plc.

ROCKET IPC REPLACEMENT VALVE INSTRUCTIONS FOR USE

Product Name: Rocket IPC Replacement Valve

Procedure Pack Product Code: R54410-00-VV

Device Image:

Device Contents:

DEVICE REF	DEVICE NAME	QTY	LEGAL MANUFACTURER DETAILS 
ZPACK005	Replacement Valve	1	ROCKET MEDICAL PLC Sedling Road, Washington, NE38 9BZ, UK www.rocketmedical.com  1639
	IPC Catheter Clamp	2	
	Scissors	1	
	Valve Cap	1	
	Slit Foam Dressing	1	
	Clear Dressing	1	
	Gauze Swab	5	
	Disposable Gloves	1	
ZBAGS137	Tiger Clinical Waste Bag	1	Abbey Polythene LTD , Unit 2, Nations Farm Curdrige Lane, Curdrige, Southampton, SO32 2BH
ZSWAB014	Skin Cleansing Alcohol Preps	4	Covidien , 15 Hampshire St. Mansfield, MA 02048

Intended Use: The Rocket IPC (Indwelling Percutaneous Catheter) Replacement Valve is intended to connect to the IPC in order to replace a damaged valve if necessary so that the drainage can continue.

Indications: For replacement of damaged IPC valve

Contraindications: Not to be used with anything else but 16fg silicone IPC valve

This valve is designed to work with silicone catheters with an internal diameter of 2.7mm (8FG), suitability to be used with the chosen catheter is at the discretion of the clinician responsible for the patients care.

WARNING: If catheter becomes damaged, clamp using supplied slide clamp or catheter clamp near to the body and seek advice with regards to preventing further damage to the catheter.

Benefit Risk: Rocket Medical Plc concludes that the overall medical benefits of the Rocket Medical IPC Replacement Valve outweigh the possible risks when used according to the intended use.

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

Performance Claims and Clinical Benefit: The IPC Drainage devices are apparatus intended to be used for external drainage of pleural effusions and abdominal ascites. The Rocket IPC Replacement Valve is intended to connect to the IPC Catheter in order to replace a damaged valve if necessary so that the drainage can continue.

Undesirable Side Effects: No side effects identified

Description of the device: The Rocket R54410-00-VV Valve Replacement Kit contains replacement valve, clamp, scissors and dressing pack.

User: The intended user is a healthcare professional proficient in the undertaking of ascites drainage, working in accordance with local and national guidelines. But also, the patient themselves can operate the device once the catheter is implanted and has undertaken the training from qualified professional.

Patient Population: The intended patient demographics are adult male and female patients requiring pleural or peritoneal ascites drainage.

Environment: *This device is for use in a clean and clear home or clinical setting.*

For Single Use Only: the device may not have structural integrity to support reuse, resulting in harm to the patient and/or user. Re-use may result in harm to patients and users from cross-contamination and infection.

Sterile: the device is supplied sterile by Ethylene Oxide (EO). Do not re-sterilise; it may not be possible to achieve the required sterility assurance level upon re-sterilisation of a used device. Re-sterilisation may also compromise the structural integrity of the device, leading to device failure.

MRI Safety: MRI compatibility has not been verified.

General Information:

- This valve is designed to work with Rocket IPC Catheter and silicone catheters with an internal diameter of 2.7mm (8FG), suitability to be used with the chosen catheter is at the discretion of the clinician responsible for the patients care.

Symbols and safety signs used on the device and labelling:

Symbol or Safety Sign	Meaning	Symbol or Safety Sign	Meaning
	Manufacturer		Authorized representative in the European Community/European Union
	Caution: Federal law restricts this device to sale by or on the order of a physician.		Do not re-use
	Date of manufacture		Consult the instructions for use
	Use-by date		Caution
	Batch code		Does not contain natural rubber latex
	Catalogue number		Medical device

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

	Sterilized using ethylene oxide		Do not use if package is damaged and consult instructions for use
	Do not resterilise		Single sterile barrier system
	Single sterile barrier system with protective packaging outside		Distributor
	Importer		

Operating Instructions:

CHECK: do not use if the packaging is open or damaged or the device is damaged. Ensure you have the full contents of the kit prior to use (see 'Device Contents').

USE:

1. If a dressing is in place over the IPC catheter remove.
2. Using an aseptic technique open the valve replacement pack on to a clean and prepared surface.
3. Put on gloves.
4. Using the small blue IPC catheter clamps attach to the silicone catheter and slide shut to ensure no air or fluid can travel either way through the catheter.
5. Using the sterile scissors provided cut through the silicone catheter approx 1cm from the end of the valve
6. Discard old valve.
7. Push fit new valve, ensuring a good fit in the catheter.
8. Remove IPC catheter clamps.
9. If required proceed with drainage.
10. Redress catheter as per usual instructions.

Warning: do not service or maintain the device whilst it is in use

Disposal: This device and its accessories used with it, should be handled and disposed of in accordance with policy of the healthcare setting and with regard to all applicable regulations, including but without limitation to, those pertaining to human health & safety and care of the environment.

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

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