ROCKET OCCLUSION PLUG FOR PLEURAL VENT

ZDOCK276 Rev. 07

Language	Page
en	2-3
cs	4-5
da	6-7
nl	8-9
fr	10-12
de	13-14
el	15-17
it	18-19
pt	20-22
cn	23-24

ROCKET MEDICAL PLC
Sedling Road, Washington,
NE38 9BZ, UK

Australian Sponsor:

ROCKET MEDICAL GmbH
Am Rosengarten 48,
15566 Schöneiche.
Germany

ROCKET MEDICAL GmbH
Am Rosengarten 48,
15566 Schöneiche.
Germany

ROCKET MEDICAL LLC
50 Corporate Park Drive #890.
East Pembroke. MA. 02359.
USA

Rocket Medical Pty Ltd
Suite 6, 157 Gordon Street,
Port Macquarie.
NSW 2444. Australia

Rocket, ROCKET, and the Rocket Medical logo are registered trademarks of Rocket Medical Plc.

ROCKET OCCLUSION PLUG FOR PLEURAL VENT INSTRUCTIONS FOR USE

Product Name: Occlusion Plug for R54565 Pleural Vent

Product Code: R54565-CP

Device Image:

Intended Use: The Pleural Vent Occlusion Plug is intended for occluding the air flow path through the R54565 Rocket Pleural Vent.

Indications: The Pleural Vent Occlusion Plug is indicated for use for patients with a suspected resolved pneumothorax in adults. This device may only be used by, or under the supervision of, appropriately trained personnel in conjunction with clinical practice guidelines.

Contraindications: The Rocket Occlusion Plug is contraindicated for indicated for patients with an active pneumothorax – use may result in a tension pneumothorax that could be fatal.

WARNING:

If pneumothorax has not resolved and the Pleural Vent is to be left in position. The occlusion cap is used to check if the pneumothorax has resolved or remains active. Do not leave occlusion cap in situ. Remove immediately after confirmation of the pneumothorax status. The patient is at risk of an induced tension pneumothorax if the occlusion cap is left in situ.

Benefit Risk: Rocket Medical plc has taken all necessary steps to ensure that residual risks associated with the use of Occlusion Plug for Pleural Vent are reduced as far as possible through application of existing state of the art techniques in the design and manufacture of these medical devices to ensure safe usage. Rocket Medical plc concludes that the overall medical benefits of Occlusion Plug for Pleural Vent, outweigh the possible risks when used according to the intended use.

Performance Claims and Clinical Benefit: The Rocket Pleural Vent system is for the treatment of spontaneous, iatrogenic simple pneumothorax whilst allowing full patient mobility during treatment. The clinical benefits of this device include:

- Resolution of pneumothorax
- Reduction in Hospital Stay
- Ambulatory Care

Undesirable Side Effects:

- Increased Pneumothorax size in PSP or SSP
- Subcutaneous emphysema
- Re-Expansion Pulmonary Oedema (Non-symptomatic)
- Pain at the insertion site
- Infection
- Bleeding
- Device dislodgement
- Device blocked or kinked
- Erythema or Itch

Description of the device: A plug to prevent flow of air through the R54565 Pleural Vent.

User: The intended user is a healthcare professional proficient in the undertaking of chest drainage, working in accordance with local and national guidelines.

Patient Population: The intended patient demographics are patients requiring drainage of air from their pleural space.

Environment: The device shall be used in a sterile field.

For Single Use Only: The device may not have structural integrity to support reuse, resulting in harm to the patient and/or user. Re-use may result in harm to patients and users from cross-contamination and infection.

Sterile: The device is supplied sterile by Ethylene Oxide (EO). Do not re-sterilise; it may not be possible to achieve the required sterility assurance level upon re-sterilisation of a used device. Re-sterilisation may also compromise the structural integrity of the device, leading to device failure.

MRI Safety: MRI compatibility has not been verified.

Symbols and safety signs used on the device and labelling:

Symbol or Safety Sign	Meaning	Symbol or Safety Sign	Meaning
	Manufacturer		Keep away from sunlight
	Authorized representative in the European Community/European Union		Keep dry
	Caution: Federal law restricts this device to sale by or on the order of a physician.		Do not re-use
	Date of manufacture		Consult the instructions for use
	Use-by date		Caution
	Batch code		Does not contain natural rubber latex
	Catalogue number		Medical device
	Sterilized using ethylene oxide		Do not use if package is damaged and consult instructions for use
	Do not re-sterilise		Single sterile barrier system
	Single sterile barrier system with protective packaging outside		Importer
	Distributor		

Operating Instructions:

CHECK: do not use if the packaging is open or damaged or the device is damaged.

USE:

1. When the signal diaphragm of the Pleural Vent (R54565) remains in the down position and the patient is asymptomatic the pneumothorax may be resolved. This must be confirmed with imaging before removal of the Pleural Vent.
2. The valve system of the Pleural Vent can be occluded using the spigot located on the occlusion plug.
3. Remove the blue valve cap covering the Pleural Vent valve.
4. Place the spigot of the occlusion plug through the valve and lock into place.
5. Confirm if the pneumothorax is resolved through medical imaging.
6. If the pneumothorax is confirmed to have resolved, the Pleural Vent may be removed.
7. If airflow through the Pleural Vent is still required remove the occlusion plug immediately.
8. Cover the valve with the blue valve cap.

Disposal: This device, its accessories and the consumables used with it, should be handled and disposed of in accordance with policy of the healthcare setting and with regard to all applicable regulations, including but without limitation to, those pertaining to human health and safety and care of the environment. Failure to do so may increase the risks of infection or other microbial hazards. Take care when handling sharps, to avoid needlestick injuries. Ensure sharps are disposed of in sharps containers.

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

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