ROCKET ASPIRATION LINE FOR PLEURAL VENT

ZDOCK277 Rev. 07

Language	Page
en	2-3
cs	4-6
da	7-9
nl	10-12
fr	13-15
de	16-18
el	19-21
it	22-24
pt	25-27
cn	28-30

ROCKET MEDICAL PLC
Sedling Road, Washington,
NE38 9BZ, UK

Australian Sponsor:

ROCKET MEDICAL GmbH
Am Rosengarten 48,
15566 Schöneiche.
Germany

ROCKET MEDICAL GmbH
Am Rosengarten 48,
15566 Schöneiche.
Germany

ROCKET MEDICAL LLC
50 Corporate Park Drive #890.
East Pembroke. MA. 02359.
USA

Rocket Medical Pty Ltd
Suite 6, 157 Gordon Street,
Port Macquarie.
NSW 2444. Australia


Rocket, ROCKET, and the Rocket Medical logo are registered trademarks of Rocket Medical Plc.

ROCKET ASPIRATION LINE FOR PLEURAL VENT INSTRUCTIONS FOR USE

Product Name: Rocket Aspiration Line for Pleural Vent

Product Code: R54565-AL

Device Image:

Intended Use: The Rocket Aspiration Line is intended for attachment to the Pleural Vent (R54565) for flushing with sterile saline and the aspiration of fluid or air.

Indications: The Rocket Aspiration Line is indicated for patients with a pneumothorax requiring flushing and aspiration. This device may only be used by, or under the supervision of, appropriately trained personnel in conjunction with clinical practice guidelines.

Contraindications: None identified.

- | | |
|-----------------|--|
| WARNING: | If Aspiration Line's 3-way tap is left in the 'off' position, this will prevent the flow of air through the device and if the pneumothorax has not resolved, the patient is at risk of an induced tension pneumothorax. Do not leave Aspiration Line in situ. Remove the Aspiration Line from the Pleural Vent device as soon as possible. |
| WARNING: | The Pleural Vent has not been tested for compatibility with medicines, nor has the clinical benefit of administering medicines through the device been established. |

Benefit Risk: Rocket Medical plc has taken all necessary steps to ensure that residual risks associated with the use of Rocket Aspiration Line for R54565 Pleural Vent are reduced as far as possible through application of existing state of the art techniques in the design and manufacture of these medical devices to ensure safe usage. Rocket Medical plc concludes that the overall medical benefits of Rocket Aspiration Line for R54565 Pleural Vent, outweigh the possible risks when used according to the intended use.

Performance Claims and Clinical Benefit: The Rocket Pleural Vent system is for the treatment of spontaneous, iatrogenic simple pneumothorax whilst allowing full patient mobility during treatment. The clinical benefits of this device include:

- Resolution of pneumothorax
- Reduction in Hospital Stay
- Ambulatory Care

Undesirable Side Effects:

- | | |
|---|----------------------------|
| - Increased Pneumothorax size in PSP or SSP | - Infection |
| - Subcutaneous emphysema | - Bleeding |
| - Re-Expansion Pulmonary Oedema (Non-symptomatic) | - Device dislodgement |
| - Pain at the insertion site | - Device blocked or kinked |
| | - Erythema or Itch |

Description of the device: The R54565-AL Rocket Aspiration Line can be connected to the R54565 Pleural Vent when flushing and aspiration is required.

User: The intended user is a healthcare professional proficient in the undertaking of chest drainage, working in accordance with local and national guidelines.

Patient Population: The intended patient demographics are patients requiring drainage of air from their pleural space.

Environment: The device shall be used in a sterile field.

For Single Use Only: The device may not have structural integrity to support reuse, resulting in harm to the patient and/or user. Re-use may result in harm to patients and users from cross-contamination and infection.

Sterile: The device is supplied sterile by Ethylene Oxide (EO). Do not re-sterilise; it may not be possible to achieve the required sterility assurance level upon re-sterilisation of a used device. Re-sterilisation may also compromise the structural integrity of the device, leading to device failure.

MRI Safety: MRI compatibility has not been verified.

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

Symbols and safety signs used on the device and labelling:

Symbol or Safety Sign	Meaning	Symbol or Safety Sign	Meaning
	Manufacturer		Keep away from sunlight
	Authorized representative in the European Community/European Union		Keep dry
	Caution: Federal law restricts this device to sale by or on the order of a physician.		Do not re-use
	Date of manufacture		Consult the instructions for use
	Use-by date		Caution
	Batch code		Does not contain natural rubber latex
	Catalogue number		Medical device
	Sterilized using ethylene oxide		Do not use if package is damaged and consult instructions for use
	Do not re-sterilise		Single sterile barrier system
	Single sterile barrier system with protective packaging outside		Importer
	Distributor		

Operating Instructions:

CHECK: do not use if the packaging is open or damaged or the device is damaged.

USE:

1. Remove the blue valve cap covering the Pleural Vent valve.
2. Place the access tip of the aspiration line through the valve and lock the cap in to place.
3. For flushing of the Pleural Vent, connect a luer lock syringe primed with sterile saline (not supplied with the device), and flush as required.
4. The flushing fluid may then be aspirated. Using the connected luer lock syringe (not supplied with the device) connected to the 3-way tap and aspirate the fluid and air as required.
5. When finished with the Aspiration line, unlock and remove from the Pleural Vent as soon as possible. Cover the valve with the blue valve cap.

Disposal: This device, its accessories and the consumables used with it, should be handled and disposed of in accordance with policy of the healthcare setting and with regard to all applicable regulations, including but without limitation to, those pertaining to human health and safety and care of the environment. Failure to do so may increase the risks of infection or other microbial hazards. Take care when handling sharps, to avoid needlestick injuries. Ensure sharps are disposed of in sharps containers.

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