Product Name: Rocket IPC Dressing Pack

ZDOCK28I Rev. 07


Language	Page
en	2-4
da	5-7
de	8-10
fr	11-13
nl	14-16

1639

ROCKET MEDICAL PLC
Sedling Road, Washington,
NE38 9BZ, UK

Australian Sponsor:

ROCKET MEDICAL GmbH
Am Rosengarten 48,
15566 Schöneiche.
Germany

ROCKET MEDICAL GmbH
Am Rosengarten 48,
15566 Schöneiche.
Germany

ROCKET MEDICAL LLC
50 Corporate Park Drive #890.
East Pembroke. MA. 02359.
USA

Rocket Medical Pty Ltd
Suite 6, 157 Gordon Street,
Port Macquarie.
NSW 2444. Australia

Rocket, ROCKET, the R device and the Rocket Medical logo are registered trademarks of Rocket Medical Plc.


ROCKET IPC DRESSING PACK

INSTRUCTIONS FOR USE

Product Name: Rocket IPC Dressing Pack

Procedure Pack Product Code: R54400-00-00

Device Image:

Device Contents: Check you have the full contents of the kit prior to use.

DEVICE REF	DEVICE NAME	QTY	LEGAL MANUFACTURER DETAILS 
ZPACK008	Slit Foam Dressing	1	ROCKET MEDICAL PLC Sedling Road, Washington, England, NE38 9BZ www.rocketmedical.com
	Clear Dressing	1	
	Gauze Swabs	5	
	Disposable Gloves	2	
ZSWAB014	Alcohol Swabs	4	See individual packaging

Intended Use: The IPC Dressing Pack is intended for procedure preparation and dressing of the wound and/or the catheter by patients fitted with an Indwelling Pleural or Peritoneal Catheter (IPC). It is intended for the dressing to be changed regularly based on clinical need. A dressing will be required for as long as the IPC is in situ. Post wound healing, the dressings are intended to protect the catheter and hold it in place.

Indications: The IPC dressing pack is indicated for patients who have been fitted with an Indwelling Pleural or Peritoneal Catheter (IPC) for the active or passive drainage of pleural effusion or abdominal ascites.

Contraindications: There are no identified contraindications for the IPC Dressing Pack.

WARNING: Do not use sharp objects such as scissors near the catheter.
If you accidentally cut or pierce your catheter, use the blue slide clamp and immediately clamp off your catheter.

1. Position your fingers between your body and the cut
2. pinch the catheter closed
3. Take the small slide clamp and slip it over the catheter then push the catheter into the small end of the clamp. This will close off the catheter preventing air entering or fluid leaking.

Contact your doctor immediately.

WARNING: Contact your nurse or doctor if pain or discomfort are prolonged or severe. Pain may be an indication of infection.

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

CAUTION: Only drain as directed by your doctor. Do not change frequency or amount of fluid drained without consulting your doctor first.

Benefit Risk:

The benefits of using the IPC Dressing Pack include:

- To palliate dyspnoea due to pleural effusion
- To provide pleurodesis
- To palliate symptoms related to recurrent ascites

Undesirable Side Effects: There are no clinical residual risks that constitute undesirable side effects

Description of the device: The IPC Dressing Pack comes packaged as a sterile kit, which includes everything needed to assist in drainage of pleural effusions or abdominal ascites and proceeding wound dressing, as well as instructions on how to use the device

User: The intended user is an appropriately trained healthcare professional or a trained caregiver in a home care setting proficient in the undertaking of external drainage of pleural effusions or abdominal ascites, working in accordance with local and national guidelines.

Patient Population: The intended patient demographics are patients requiring active or passive drainage of pleural effusions or abdominal ascites. The majority of patients requiring these devices are over the age of 35 and have been diagnosed with, but not limited to, lung cancer, mesothelioma, ovarian cancer or breast cancer.

Environment: This device is for use in a clean and clear home community or clinical setting.

For Single Use Only: the device may not have structural integrity to support reuse, resulting in harm to the patient and/or user. Re-use may result in harm to patients and users from cross-contamination and infection.


Sterile: the device is supplied sterile by Ethylene Oxide (EO). Do not re-sterilise; it may not be possible to achieve the required sterility assurance level upon re-sterilisation of a used device. Re-sterilisation may also compromise the structural integrity of the device, leading to device failure.

MRI Safety: MRI compatibility has not been verified.

Symbols and safety signs used on the device and labelling:

Symbol or Safety Sign	Meaning	Symbol or Safety Sign	Meaning
	Manufacturer		Keep away from sunlight
	Authorized representative in the European Community/European Union		Keep dry
	Caution: Federal law restricts this device to sale by or on the order of a physician.		Do not re-use
	Date of manufacture		Consult the instructions for use
	Use-by date		Caution
	Batch code		Does not contain natural rubber latex
	Catalogue number		Medical device
	Sterilized using ethylene oxide		Do not use if package is damaged and consult instructions for use

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

	Do not resterilise		Single sterile barrier system
	Distributor		Single sterile barrier system with protective packaging outside
	Importer		

Operating Instructions:


CHECK: do not use if the packaging is open or damaged or the device is damaged. Ensure you have the full contents of the kit prior to use (see 'Device Contents')

USE:


1. Use a clean, clear workspace, ideally a table or bench top.

Important: Wash your hands with warm soapy water for at least 1 minute before you begin.

2. Remove the dressings from around the catheter site, if you notice anything unusual such as pain, redness, swelling or free fluid, finish this procedure first and then contact your doctor or nurse immediately.
3. Open your procedure pack and position the self-adhesive dressing to one side. Position the blue wrap on your clean workspace with the tape facing upwards. Carefully remove tape and unfold the blue wrapping making sure the contents remain in the middle of the sterile wrapping.
4. Put the gloves on. It is important that these remain sterile. Pick the first one up by the inside of wrist fold and place on the opposite hand, repeat this for the other hand. Your bare hand should not come into contact with the outer surface of the glove.
5. Clean around the catheter exit site with an alcohol pad (**Fig.1**).
6. Place the foam catheter pad around the catheter (**Fig.2**).
7. Wind the catheter carefully into loops and place over the foam (**Fig.3**). Cover the catheter with the gauze pads (**Fig.4**).

8. Take the wound dressing and peel away the paper backing from one side (**Fig.5**).
9. Position the dressing centrally over the gauze pads and carefully press down (**Fig.6**), then peel the other half of the dressing away pressing firmly all the time. Hold one of the coloured tabs on the dressing and pull the cover sheet away from the rest of the dressing, smooth and press until secure (**Fig.7**).


Disposal: This device, its accessories and the consumables used with it, should be handled and disposed of in accordance with policy of the healthcare setting and with regard to all applicable regulations, including but without limitation to, those pertaining to human health and safety and care of the environment. Failure to do so may increase the risks of infection or other microbial hazards. Take care when handling sharps, to avoid needlestick injuries. Ensure sharps are disposed of in sharps containers.

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

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