**Product Names:** Oocyte Aspiration Pump Filter Sets

## ZDOCKI94 Rev. 19

| Language | Page  |
|----------|-------|
| EN       | 1 - 4 |
| CN       | 5 - 7 |
| CZ       | 8-11  |
| DA       | 12-15 |
| DE       | 16-19 |
| EL       | 20-23 |
| ES       | 24-27 |
| FI       | 28-31 |
| FR       | 32-35 |
| IT       | 36-39 |
| NL       | 40-43 |
| NO       | 44-47 |
| PL       | 48-51 |
| SV       | 52-55 |
| VI       | 56-59 |

**1639**

**ROCKET  
MEDICAL PLC**  
Sedling Road,  
Washington,  
United Kingdom  
NE38 9BZ

**ROCKET MEDICAL GmbH**  
Am Rosengarten 48,  
15566 Schöneiche.  
Germany

**ROCKET MEDICAL GmbH**  
Am Rosengarten 48,  
15566 Schöneiche.  
Germany

**ROCKET MEDICAL LLC**  
50 Corporate Park Drive #890.  
East Pembroke. MA. 02359.  
USA

Australian Sponsor:

**Rocket Medical Pty Ltd**  
Suite 6, 157 Gordon Street,  
Port Macquarie.  
NSW 2444. Australia

The R device and the Rocket Medical logo are registered trademarks of Rocket Medical Plc.

### Product Codes:

R57686 Filter Set for Oocyte Aspiration Pump

R57687 Filter Set for Cook K-MAR-5200 Pump

### Device Image:

Figure 1: Filter Set for Oocyte Aspiration Pump

### Intended Use:

**R57686 Filter Set for Oocyte Aspiration Pump:** Intended for the connection of oocyte aspiration needles to the Rocket Medical Oocyte Aspiration Pump.

**R57687 Filter Set for Cook K-MAR-5200 Pump:** Intended for the connection of oocyte aspiration needles to the Cook K-MAR-5200 Pump.

**Indications:** Indicated for use as part of the treatment for infertility relating to IVF and other related gynaecological procedures. This device may only be used by, or under the supervision of, appropriately trained personnel in conjunction with clinical practice guidelines.

**Contraindications:** Do not allow to become contaminated with fluid. If the device becomes contaminated replace immediately.

**WARNING:** The hydrophobic filter set is specially designed to prevent ingress of fluid to the pump. If the filter material becomes contaminated with fluid it will occlude thereby protecting the pump mechanism. Immediately replace the filter set if it becomes contaminated and between every patient.

**WARNING:** Do not allow to become contaminated with fluid. If the device becomes contaminated replace immediately

**WARNING:** Do not service or maintain the device whilst it is in use

**Precaution:** **Potential for DEHP Exposure from use of this device.** DEHP may leach from this device, which may cause reproductive harm in male neonates, pregnant women carrying male fetuses, and peripubertal males. The extent of exposure largely depends upon the medical treatments administered and the duration of the treatment. Lack of research in humans means it is difficult to predict the adverse effects of DEHP because certain animal models may not apply to humans. Most studies utilise mice and rats, however breakdown of DEHP in the human body differs to the mechanism observed in these animals. These effects are also seen only at levels far in excess of normal human exposure. (For further details please see 'POL58' at [www.rocketmedical.com](http://www.rocketmedical.com).)

**Benefit Risk:** Rocket Medical plc has taken all necessary steps to ensure that residual risks associated with the use of Rocket Medical filter sets are reduced as far as possible through application of existing state of the art techniques in the design and manufacture of these medical devices to ensure safe usage. Rocket Medical plc concludes that the overall medical benefits of Rocket filter sets, outweigh the possible risks when used according to the intended use.

*In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at [pncf@rocketmedical.com](mailto:pncf@rocketmedical.com) as well as to the competent health authority in the country that the user/patient resides.*

**Performance Claims/Characteristics and Clinical Benefit:** The performance claim of the filter sets is to be able to protect the oocyte aspiration pump from contamination. The clinical benefit is treatment for infertility relating to IVF and other related gynaecological procedures.

**Undesirable Indirect Side Effects:** bleeding, infection, pain and endometrioma perforation/ injury of pelvic structures

**Description of the device:** This device is a 2.5m lightweight patient tube set with an integral Luer lock connector to ensure a secure, reliable connection to all standard Luer lock fitting aspiration needle sets. This filter set is supplied sterile and is intended for single-use only.

**User:** The intended user is a healthcare professional proficient in the undertaking of oocyte aspiration, working in accordance with local and national guidelines. The device may also be used by a trainee under the supervision of a healthcare professional proficient in the undertaking oocyte aspiration.

**Patient Population:** The intended patient demographics are adult patients undergoing oocyte aspiration for the treatment of infertility with IVF or other gynaecological patients such as those undergoing oocyte retrieval for cryopreservation or those with ovarian hyperstimulation syndrome.

**Environment:** For use in a clinical setting.

**For Single Use Only:** the device may not have structural integrity to support reuse, resulting in harm to the patient and/or user. Re-use may result in harm to patients and users from cross-contamination and infection.

**Sterile:** the devices are supplied sterile by Ethylene Oxide (EO). Do not re-sterilise; it may not be possible to achieve the required sterility assurance level upon re-sterilisation of a used device. Re-sterilisation may also compromise the structural integrity of the device, leading to device failure.

**MRI Safety:** MRI compatibility has not been verified.

**Symbols and safety signs used on the device and labelling:**

| Symbol or Safety Sign                                                               | Meaning                                                                               | Symbol or Safety Sign                                                               | Meaning                                                           |
|-------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------|
|  | Manufacturer                                                                          |  | Keep away from sunlight                                           |
|  | Caution: Federal law restricts this device to sale by or on the order of a physician. |  | Keep dry                                                          |
|  | Date of manufacture                                                                   |  | Do not re-use                                                     |
|  | Use-by date                                                                           |  | Consult the instructions for use                                  |
|  | Batch code                                                                            |  | Caution                                                           |
|  | Catalogue number                                                                      |  | Does not contain natural rubber latex                             |
|  | Sterilized using ethylene oxide                                                       |  | Medical device                                                    |
|  | Do not re-sterilise                                                                   |  | Do not use if package is damaged and consult instructions for use |

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at [pncf@rocketmedical.com](mailto:pncf@rocketmedical.com) as well as to the competent health authority in the country that the user/patient resides.

|                                                                                   |                                                                 |                                                                                   |                               |
|-----------------------------------------------------------------------------------|-----------------------------------------------------------------|-----------------------------------------------------------------------------------|-------------------------------|
|  | Single sterile barrier system with protective packaging outside |  | Single sterile barrier system |
|  | Distributor                                                     |  | Importer                      |

### Operating Instructions:

**CHECK:** do not use if the packaging is open or damaged or the device is damaged.

**USE:**

1. Carefully follow the instructions for assembling and setting up the Oocyte Aspiration Pump before attempting to attach the Patient Connection Set.

**For Rocket Medical Oocyte Aspiration Pump follow instructions 2-3 below regarding powering on the device. For Cook K-MAR-5200 follow the powering on instructions provided in your user manual.**

2. Connect the correct mains lead to the 12V power supply module and connect to a 110-240VAC 50-60Hz electrical supply. The power supply module is self-regulating to generate 12V DC.
3. Turn the power switch 0-I (rear panel) to on. The green LED (indicated) on the front panel will illuminate and the main display will display the pre-set value of -100mmHg.

### Attaching the filter set

4. Using aseptic technique, un-pack the Filter Set and attach the short tubing length, proximal to the filter, to the pump spigot (indicated right for Rocket Medical Oocyte Aspiration Pump).
5. Pass the longer, patient connection tube into the operative field. The filter set is for single patient use and must be replaced for each patient.

**Set the desired vacuum level. For Rocket Medical Oocyte Aspiration Pump follow points 6-9 below. For Cook K-MAR-5200 follow the instructions provided in your user manual.**

6. On Start-Up, the pump is pre-set to produce -100mmHg vacuum. To increase this value, use a circular motion to lightly draw a finger clockwise over the control dial in the direction of the arrows. Refer to Table 1 for recommend values for various needle and tube set combinations.
7. The dial will sense the motion and increase the display in -1mmHg increments.
8. To decrease the set value, use an anticlockwise circular motion.
9. When the desired value is reached, hold the finger still for 2 seconds and an audible 'beep' will confirm the new setting. If the new setting is not correctly saved, the value will return to its previous set value in 10 seconds. Note that if the finger is held still for too long (8-10 seconds), 3 beeps will sound and the pump will enter service mode, causing the displayed value to fluctuate during operation. Service mode should only be used by trained engineers during pump repair or service, and not during clinical use. If service mode has been activated in error, hold a finger still on the outer zone of the control dial for another 8-10 seconds to return to the normal display mode, Two beeps will occur to confirm the change.

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at [pncf@rocketmedical.com](mailto:pncf@rocketmedical.com) as well as to the competent health authority in the country that the user/patient resides.

|          | Single Lumen             |     | Double Lumen |     |
|----------|--------------------------|-----|--------------|-----|
|          | 16G                      | 17G | 16G          | 17G |
| Tube Set | Recommended Vacuum -mmHg |     |              |     |
| 55cm     | 80                       | 110 | 130          | 150 |
| 70cm     | 90                       | 130 | 150          | 170 |
| 90cm     | 100                      | 150 | 170          | 190 |

Table 1: Recommended Vacuum Settings for Needle Sizes and Tubing Length.

**Disposal:** This device, its accessories and the consumables used with it, should be handled and disposed of in accordance with policy of the healthcare setting and with regard to all applicable regulations, including but without limitation to, those pertaining to human health and safety and care of the environment. Failure to do so may increase the risks of infection or other microbial hazards.

*In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at [pncf@rocketmedical.com](mailto:pncf@rocketmedical.com) as well as to the competent health authority in the country that the user/patient resides.*

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