Product Name: FBS Sampling Wand**ZDOCKI85 Rev. 14**

Language	Page
en _____	2-6
nl _____	7-11
fr _____	12-16
da _____	17-21
es _____	22-26
de _____	27-31
fi _____	32-36
sv _____	37-41
no _____	42-46

ROCKET MEDICAL PLC

Sedling Road, WASHINGTON,
United Kingdom, NE38 9BZ

ROCKET MEDICAL GmbH

Am Rosengarten 48,
15566 Schöneiche.
Germany

**ROCKET MEDICAL LLC**

50 Corporate Park Drive #890.
East Pembroke. MA. 02359.
USA

Rocket Medical Pty Ltd

Suite 6, 157 Gordon Street,
Port Macquarie.
NSW 2444. Australia

Australian Sponsor:

The R device and the Rocket Medical logo are registered trademarks of Rocket Medical Plc. Amnilume™ is a trademark of Rocket Medical Plc.

Product Name: Rocket Medical Fetal Blood Sampling (FBS) Sampling Wands with Capillary Tubes

Product Codes: R57021-55-SW, R57021-85-SW and R57021-10-SW

Device image

Device Contents: Check you have the full contents of the kit prior to use.

DEVICE REF	DEVICE NAME	QTY	LEGAL MANUFACTURER DETAILS 
R57021-55-SW	FBS Sampling Wand	1	ROCKET MEDICAL PLC , Sedling Road, Washington, United Kingdom, NE38 9BZ www.rocketmedical.com
	55µL Capillary Tube	5	RADIOMETER MEDICAL APS , Åkandevvej 21, 2700 København, Denmark radiometer@radiometer.dk
R57021-85-SW	FBS Sampling Wand	1	ROCKET MEDICAL PLC , Sedling Road, Washington, United Kingdom, NE38 9BZ www.rocketmedical.com
	85µL Capillary Tube	5	RADIOMETER MEDICAL APS , Åkandevvej 21, 2700 København, Denmark radiometer@radiometer.dk
R57021-10-SW	FBS Sampling Wand	1	ROCKET MEDICAL PLC , Sedling Road, Washington, United Kingdom, NE38 9BZ www.rocketmedical.com
	100µL Capillary Tube	5	SIEMENS HEALTHCARE DIAGNOSTICS INC , 511 Benedict Avenue, Tarrytown, New York 10591-5097 USA

Intended Use: The sampling wand is intended to be used to make a small incision in the fetal scalp and to hold a co-supplied, dimensionally matched capillary tube, in order to collect fetal blood sample(s) for subsequent analysis.

Indications: The fetal blood sampling wands are indicated for use for patients requiring fetal blood analysis (for pH/blood gas) during labour. Dilation of at least 3cm is required, with ruptured membranes.

Contraindications:

Fetal blood sampling is contraindicated:

- when there is acute fetal compromise (for example, fetal bradycardia of >3 minutes)
- when there is an acute event such as cord prolapse, uterine rupture or suspected placental abruption
- where the full clinical picture suggests that the birth should be expedited, for example, in cases of maternal sepsis or unstable pre-eclampsia

as this would delay an urgent requirement to deliver the fetus, leading to deterioration in the fetal condition.

FBS is also contraindicated:

- in cases of placenta previa as this may cause serious bleeding.
- in maternal thrombocytopenia as this may lead to increased bleeding, bruising, or abnormal blood clotting.
- when the fetus is premature (less than 34+0 weeks) as this may lead to or worsen existing fetal anaemia.
- where the mother is a confirmed carrier of haemophilia and the fetus is either affected or of unknown status as this may lead to increased bleeding.
- where there is a risk of maternal-fetal infection transfer (for example, active herpes simplex virus, HIV and hepatitis C). Use caution and refer to your local clinical policies and procedures in cases when the woman is positive for Group B streptococcus, gonorrhoea, hepatitis A, hepatitis B or syphilis.
- in cases of facial presentation as there may be permanent facial scarring or injury to the eye.
- when identification of the presenting part is uncertain as an incision may be made into the fetal genitals, face or fontanelles, leading to permanent injury to the eye, loss of CSF or scarring.

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

WARNING: Ensure close contact between amnioscope and fetal scalp is maintained during the sampling procedure to prevent leakage of amniotic fluid into the blood sample, as this could lead to a falsely reassuring result.

WARNING: Do not use this kit without prior verification that it is suitable for use with the analyser.

WARNING: There is uncertainty about the usefulness of FBS in women with fever, sepsis or suspected sepsis, since FBS may be falsely reassuring. In these cases, the decision to perform FBS or results obtained from FBS should be discussed with a consultant obstetrician.

WARNING: This device contains 2 lithium-based primary cells.
Do NOT attempt to disassemble the device.
Do NOT attempt to recharge or short circuit the batteries.

CAUTION: Avoid fetal face, fontanelles and genitalia. Rocket Medical FBS blades have a maximum depth of penetration of 2.1mm

CAUTION: Do not use Ethyl Chloride Spray – this may adversely affect the internal surface of the plastic amnioscope.

CAUTION: Ensure that both the obturator and Amnilume™ are in position in the amnioscope before inserting into the patient. If the Amnilume™ is not in place, the obturator can protrude beyond the end of the amnioscope, which may cause patient tissue to become trapped between the amnioscope and the obturator.

CAUTION: Carefully follow the analyser manufacturer's instructions to determine the capillary tube required. Failure to use correct blood volume, select correct analysis mode or the use of clotted or contaminated samples may lead to errors in results or rejection of the whole sample.

CAUTION: Do NOT use the sample if the blood clots.

CAUTION: Avoid distributing air through the sample; samples containing air bubbles may be rejected by the analyser.

Benefit Risk: Rocket Medical plc has taken all necessary steps to ensure that residual risks associated with the use of Rocket Medical FBS are reduced as far as possible through application of existing state of the art techniques in the design and manufacture of these medical devices to ensure safe usage. Rocket Medical plc concludes that the overall medical benefits of Rocket Medical FBS outweigh the possible risks when used according to the intended use.

Performance Claims/Characteristics and Clinical Benefit:

The intended performance claims for the device are:

- To make a small incision in the fetal scalp
- Device is safe to use
- To collect blood suitable for subsequent analysis

The intended clinical benefit is:

- The determination of fetal blood pH level to aid in decision making during labour

Undesirable Side Effects: Fetal side effects include: bleeding, haematoma, infection at incision point.

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

Description of the Device: The sampling wand is a device that combines the means to puncture the fetal scalp and to hold the capillary tube for blood sample collection by capillary action. In preparation for sample collection, a capillary tube is placed into the holder of the Sampling Wand.

Following scalp preparation, the Sampling Wand is then inserted into the vagina through an amnioscope (not provided with FBS Sampling Wands, but available separately). When the fetal scalp has been punctured and a droplet/pool of fetal blood has formed, the capillary tube is presented to the droplet of blood, which enters the tube by capillary action. The capillary tube is inserted into a blood gas analyser, which will provide the pH/blood gas measurement. The measurement is interpreted using clinical guidelines.

One procedure pack may take up to five samples through presentation of multiple capillary tubes to the blood droplet.

The blade is not intended to be reused since it may have become blunt and reuse could also present a risk of infection.

User: the intended user is a healthcare professional proficient in the undertaking of fetal blood sampling, working in accordance with local and national guidelines. The device may also be used by a trainee under the supervision of a healthcare professional proficient in the undertaking of FBS.

Patient Population: the intended demographic for this device is pregnant patients during labour.

Environment: to be used in a clinical setting.

For Single Use Only: the device may not have structural integrity to support reuse, resulting in harm to the patient and/or user. Re-use may result in harm to patients and users from cross-contamination and infection.

Sterile: the devices are supplied sterile by Ethylene Oxide (EO). Do not re-sterilise; it may not be possible to achieve the required sterility assurance level upon re-sterilisation of a used device. Re-sterilisation may also compromise the structural integrity of the device, leading to device failure.

MRI Safety: This device is not intended to be routinely found in an MRI environment and compatibility has not been verified.

General Information:

Rocket Medical supplies a number of FBS Sampling Wand kits with different capillary tubes. Please check your analyser user manual to confirm capillary tube compatibility. Where in doubt, please contact your analyser manufacturer representative. Use the lowest volume compatible tube size to minimise the volume of blood required for collection.

Product Code	Tube Volume / μ l	Tube Length	Tube OD	Tube ID	Heparin Type	Concentration /IU/ml
R57021-55-SW*	55	59mm +2.0 - 0.0mm	2.05 $\pm$ 0.1mm	1.25 $\pm$ 0.05mm	Na	240
R57021-85-SW*	85	87 mm +/- 1 mm	2.10 mm + 0.05 mm / - 0.10 mm	1.25 mm + 0.10 mm / - 0.05 mm	Na	240
R57021-10-SW	100	85.0mm	2.3 $\pm$ 0.05mm	1.20 $\pm$ 0.05mm	Ca/Na/Li	200

*Not suitable for lactate measurement on Radiometer analysers as this measurement requires a balanced heparin tube.

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

Symbols and safety signs used on the device and labelling:

Symbol or Safety Sign	Meaning	Symbol or Safety Sign	Meaning
	Manufacturer		Keep away from sunlight
	Authorized representative in the European Community/European Union		Keep dry
	Caution: Federal law restricts this device to sale by or on the order of a physician.		Do not re-use
	Date of manufacture		Consult the instructions for use
	Use-by date		Caution
	Batch code		Does not contain natural rubber latex
	Catalogue number		Medical device
	Sterilized using ethylene oxide		Do not use if package is damaged and consult instructions for use
	Do not resterilise		Single sterile barrier system
	Single sterile barrier system with protective packaging outside		Importer
	Distributor		

Operating Instructions:


CHECK: do not use if the packaging is open or damaged or the device is damaged. Ensure you have the following consumables and accessories available to perform the fetal blood sample (the following are not provided with FBS Sampling Wands):

- Amnioscope and light source
 - Rocket Medical recommends the use of a Rocket Medical amnioscope fitted with Amnilume™ Single Use Light Source (R57016-LUME-NS or R57016-PLUS-NS).
- Gauze swabs
- Sponge holding forceps
- Swab stick
- Sterile paraffin/petroleum jelly
- Capillary tube holder (optional)
- Capillary tube caps

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

USE:

1. Using aseptic technique, remove the device from its outer bag.
2. Select a capillary tube from the pack and load into the holder. Locate the end of the tube into the holder at (1) and then advance into the rear tube grip at (2). Ensure the tube is secure before sampling.
3. Follow local hospital procedure to prepare the perineum and vagina and complete a vaginal examination to establish the cervical dilatation and station of the fetal head.
4. Insert amnioscope into the vagina and position light source as per instructions for use of the chosen device.
5. Position the amnioscope to obtain the best possible seal against the fetal presentation scalp to prevent leakage of amniotic fluid into the blood sample.
6. Clean the puncture area using gauze swabs and sponge holding forceps.
7. With a clean swab stick, apply a *thin* layer of white soft paraffin/petroleum jelly to the puncture area to promote capillary blood aggregation. **Remove excess with a new, dry swab stick before performing the skin puncture.**
8. The sampling wand comes pre-fitted with an FBS blade. Lightly holding the shaft and rear finger grip (A), push the grip forward to expose the blade (B). **Ensure that the blade is locked forward before use.**
9. Perform the skin puncture using either technique (i) make a vertical stab, rotate the blade 90° and make a second stab to form a cross or (ii) make a vertical stab, then angle the blade upwards and extend using the front edge of the blade to produce a single 4mm incision. **Always wait for a good sized (3- 4mm) bead to form on the surface before collection.** Never slash at the scalp and avoid multiple small stabs.
10. To collect the blood sample, depress the button on the wand shaft to retract the FBS blade, and angle the capillary tube into the bead.
11. Rapidly collect the blood sample, filling the tube completely, **without air bubbles**. Collection needs to be expedited to avoid blood clotting.
12. Twist lightly to remove the capillary tube from sampling wand and cap both ends. Take to the analyser **immediately**. Follow the analyser manufacturer's instructions as regards sample mixing.
13. Always ensure that the **clean**, non-paraffin/petroleum jelly exposed end is presented to the analyser. Note: some analysers require a 'clot catcher' (not provided). Please refer to the manufacturer's instructions for guidance on when a clot catcher is needed.
14. Always ensure that the **correct** measuring mode is selected before performing the analysis.

Post-sampling:

15. Haemostasis is accomplished by applying pressure to puncture site with a swab for one full contraction.

Disposal: This device and the consumables used with it, should be handled and disposed of in accordance with policy of the healthcare setting and with regard to all applicable regulations, including but without limitation to, those pertaining to human health and safety and care of the environment. Failure to do so may increase the risks of infection or other microbial hazards. Take care when handling sharps, to avoid needlestick injuries. Ensure sharps are disposed of in sharps containers.

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

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